

OPERATION

DAILY:

- Run pool equipment and oxygen generator a minimum of 12 hours daily, during daylight hours.
- Running the system longer will provide even better quality water. Running the system continuously will provide the highest quality water.

WEEKLY:

- Test pool water with the Oxygen Pools 4 in 1 Test Kit to verify levels.
- Add recommended dose of Formula "O" as directed on pail

MONTHLY:

- Bring water sample to dealer to test for:
 - pH
 - Total Alkalinity
 - Calcium Hardness
 - Copper

Treat water according to test results.

- Clean filter by backwashing or washing cartridge

ANNUALLY:

- Check clear vinyl tubing and clamps on manifold to see if replacement is needed due to discoloration and degradation.
- Check inside Dynamic Oxygen Generator to see if vinyl tubing and check valves need replacement due to discoloration and degradation

WINTERIZATION

Flex Hose:

- Plug skimmer and return. Drain Oxygen Generator at the same time you drain the rest of your pool equipment.
- Remove entire manifold with generator by unclamping the flex hose from one end and unthreading the grey fitting out of the filter. Store in warm, dry place.
- Cover opening to filter.

PVC Pipe:

- Drain Oxygen Generator at the same time you drain the rest of your pool equipment.
- Remove the Oxygen Generator manifold from the manifold by loosening the unions on each end.
- Tape the openings over open ends of the elbows.
- Store the generator, and manifold in a warm, dry place.